



7. Turn the laser on
8. Align the optical components starting with the component nearest to the laser. When it is aligned, move the first beam stopper behind the third optical component. Repeat this procedure until the entire optical system is aligned. It is important that the laser beams be limited to one new component at a time until the system is aligned. This will minimize uncontrolled reflection during the alignment procedure

### **IMPORTANT:**

**DO NOT REMOVE YOUR SAFETY EYEWEAR DURING THE ALIGNMENT PHASE. IF YOU CANNOT SEE A FAINT IMAGE OF THE BEAM YOU HAVE THE WRONG OPTICAL DENSITY EYEWEAR. TURN OFF THE LASER AND OBTAIN THE CORRECT OPTICAL DENSITY**

### **Laser work**

1. Change the alignment goggles with normal operating goggles
2. Ensure everyone in the room wears that appropriate laser safety eyewear
3. Increase beam power (if necessary) and complete the assigned task
4. Always use the lowest beam power necessary for the procedure

### **Shut Down Procedure**

1. Turn off the laser
2. Remove your laser safety eyewear and place it in the proper storage area
3. Allow the laser to cool down and turn off the cooling water
4. Remove the key from the laser interlock switch
5. Turn off the laboratory interlock system and the laser light
6. Return the key to the laser supervisor/proper location

### **In case of emergency**

1. If possible shut down the laser by using the emergency button or by removing the laser key
2. If shut down of the laser is not possible alert everyone to leave the lab and leave the lab yourself
3. In case of fire follow the emergency procedure established in your lab
  - Familiarize yourself with the location of fire extinguishers in your area
  - Familiarize yourself with fire exits in your area
  - Attempt to extinguish the fire only if you can do it safely using one fire extinguisher
  - If not successful leave the room using the pre-established exit
  - Set the fire alarm
  - Inform U of T fire department of any fire doesn't matter how small